

EZ-PEAZY™

CPAPology®

Powerful Enzymatic Cleaning. Just Add Water.

Multi-Enzyme Cleaning Action

Specially formulated enzymes break down and lift dried organic buildup, including mucus, skin oils and protein residue.

Highly Concentrated 1:127 Formula

A little goes a long way – just 15 mL (½ oz.) makes 1.9 L (2 qt.) of ready-to-use cleaning solution.

Low-Foaming & Easy to Rinse

Low-foaming action allows for thorough rinsing and makes it ideal for cleaning CPAP masks and accessories.

Gentle on CPAP Equipment

Non-corrosive and phosphate-free, with a non-caustic formula designed for regular CPAP cleaning.

Versatile Cleaning


Suitable for a wide range of CPAP masks, accessories and cleaning devices.



Dilution 1:127 (15 mL per 1.9 L of water)



1 Organic Buildup
Everyday use can leave behind oils, mucus and protein residue.



2 Enzymes Break It Down
Specially formulated enzymes target and break down organic residue into smaller particles.



3 Residue Rinses Away
Loosened particles rinse away easily, leaving your CPAP equipment clean and fresh.

Product	Part #	UPC #
EZ-PEAZY™ Enzymatic CPAP Soap - 1 L (ea.)	CY-K8630	663427408202
EZ-PEAZY™ Enzymatic CPAP Soap - 100 mL (Case of 24)	CY-K8632	663427450461

KEGO

1-866-862-7328 | www.kegocorp.com | www.kegousa.com

www.cpapology.ca